



**Molecular
Dimensions**

A CALIBRE SCIENTIFIC COMPANY

Data Sheet

The ANGSTROM Additive Screen
Optimisation, Additive Screen

The ANGSTROM Additive Screen HT-96 and FX-96 Pre-Filled Plate MD1-100 and MD1-100-FX

Key Info

The ANGSTROM Additive Screen is a 96-condition additive screen incorporating 30 different polyols to aid protein crystallization optimization and flash-cooling of protein crystals. MD1-100 is presented as 96 x 1 mL conditions and MD1-100-FX is presented as 96 x 100 μ L.

- For easy optimization of crystal conditions.
- Contains 30 different polyols at three different concentrations.
- Cryoprotect and enhance the stability of your protein and improve crystal quality.
- Suitable for both hand and robotic dispensing.

Introduction

The ANGSTROM Additive Screen provides an easy way to optimize your crystallization condition. Once initial hits have been obtained, test out an extensive array of 30 polyols facilitated to improve the quality of your crystals and give you optimum control over your optimization experiment.

The traditional method of cryoprotection by soaking a crystal in cryoprotectant can result in damaged or lost crystals; The ANGSTROM Additive Screen contains derivatives of glycols (Figure 1, a), carbohydrates and PEGs that will be integrated to the condition early to help with cryoprotection.

The concentrations of cryoprotectant required to provide adequate cryoprotection have been marked on the Table of screen contents (page 3). Some cryoprotectants found in the ANGSTROM screen are as potent as glycerol (Figure 1).




Figure 1. (a) Glycol derivatives. From top to bottom: 1,2,3-propanetriol (glycerol), 1,2,4-butanetriol, 1,2,6-hexanetriol, 1,5-pentanediol and 1,1,1-tris(hydroxymethyl)propane. In addition to being well-suited as crystallization reagents, these five polyols are cryo-protectants when used at concentrations as low as 20–25% (w/v). (b) Crystals of endosomal sorting complex required for transport (ESCRT)-I. Ten percent of the ANGSTROM screen was added to the reservoirs of a crystallization plate pre-filled with 96 repeats of the initial condition. Several hits were observed including the one shown here with glucose as additive (final conc. 3%, w/v). Unpublished results obtained during the early stage of the ANGSTROM screen development. Work of Nicolas Soler (LMB).





**Molecular
Dimensions**

A CALIBRE SCIENTIFIC COMPANY

Data Sheet

The ANGSTROM Additive Screen
Optimisation, Additive Screen

The ANGSTROM Additive Screen HT-96 and FX-96 Pre-Filled Plate MD1-100 and MD1-100-FX

Tips for Use

The ANGSTROM Additive Screen is stored frozen at -20°C. Thaw additive screen prior to use and spin down ~1000G. Use the ANGSTROM Additive Screen from either the deep-well block or the Costar (FX) plates. The deep-well blocks contain approx. 1mL of additive screen and the Costar FX plates contain 100 µL of additive screen. We recommend aliquoting from the deep-well block into smaller volumes or directly into the reservoirs of MRC plates for example. This is to avoid repeated freeze-thaw cycles and contamination (especially if working in a multi- user lab).

Use of the ANGSTROM additive MD1-100 Deep-well block screen.

Dispense directly from the deep-well block into a storage plate to use on your robot, or directly from the block if using by hand.

Add The ANGSTROM Additive Screen into the reservoir of a crystallization plate.

If starting from a crystallization plate (e.g. MRC plate) it is recommended that you pre-dispense 8 µL of The Angstrom Additive Screen into the reservoir and then dispense 72 µL of your crystallization screen on top. Place plate on orbital shaker to mix prior to setting protein. Then continue with dispensing your protein etc.

Add The ANGSTROM Additive Screen directly into the drop.

Dispense directly into the sample well of a crystallization plate, for example, 500 nL of protein + 500 nL crystallization condition + 100 nL of the additive screen. You can vary the ratios of protein drop to crystallization drop and additive drop.

Formulation Notes

The Angstrom Additive Screen reagents are formulated using ultrapure water (>18.0 MΩ) and are sterile filtered using 0.22 µm filters. No preservatives are added. Final pH may vary from that specified on the datasheet. Molecular Dimensions will be happy to discuss the precise formulation of individual reagents. Individual reagents and stock solutions for optimization are available from Molecular Dimensions.

Unit 5A, R-Evolution @
The Advanced Manufacturing Park
Selden Way, Rotherham, S60 5XA



**Molecular
Dimensions**

Enquiries regarding The ANGSTROM Additive Screen formulation, interpretation of results or optimization strategies are welcome.

Please e-mail or phone your query to Molecular Dimensions.

Contact and product details can be found at:

www.moleculardimensions.com

Manufacturer's safety data sheets are available from our website.

The conditions shown on this datasheet may differ from those shown on previous versions of the datasheet due to the discontinuation for raw material supply for the following:
Pentaerythritol ethoxylate (3/4 EO/OH)

References

1. Gorrec, F. Protein crystallization screens developed at the MRC Laboratory of Molecular Biology, Drug Discov Today (2016), <http://dx.doi.org/10.1016/j.drudis.2016.03.008>

The ANGSTROM Additive Screen has been designed and developed by Fabrice GORREC, in collaboration with the scientists at the Medical Research Council Laboratory of Molecular Biology (LMB) at Cambridge and is manufactured exclusively under license by Molecular Dimensions Limited.

Limited Use and Restrictions: Products sold by Molecular Dimensions Ltd. or its affiliates or authorized distributors and information relating to same are intended for research use only in crystal growth and optimization of crystal growth following use of the product by the purchaser and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, including resale or use in manufacture. The license to use The Angstrom Additive Screen specifically excludes any rights to use the product information for the manufacture of the product or derivatives thereof, or distribute, transfer, or otherwise provide access to such information to any third party for any purpose or use.

Ordering Information

PN:	Description:	Pack Size:
MD1-100	ANGSTROM HT-96	96 x 1 mL
MD1-100-FX	ANGSTROM FX-96 pre-filled plate	96 x 100 µL
Single Reagents		
MDSR-100-well number	ANGSTROM HT-96 single reagent	100 mL

enquiries@moleculardimensions.com
+44 (0)114 224 2257

calibre
SCIENTIFIC



**Molecular
Dimensions**

A CALIBRE SCIENTIFIC COMPANY

Data Sheet

The ANGSTROM Additive Screen
Optimisation, Additive Screen

The ANGSTROM Additive Screen HT-96 and FX-96 Pre-Filled Plate HT-96 MD1-100 | MD1-100 FX | Conditions A1 to H12

Well #	Conc.	Units	Reagent	Cryoprotecting	Well #	Conc.	Units	Reagent	Cryoprotecting
A1	100	%	Ultrapure water	-	E1	15	% w/v	D-Glucose	-
A2	100	%	Ultrapure water	-	E2	30	% w/v	D-Glucose	-
A3	100	%	Ultrapure water	-	E3	60	% w/v	D-Glucose	Yes
A4	0.5	M	Lithium sulphate	-	E4	2.5	% w/v	<i>myo</i> -Inositol	-
A5	1	M	Lithium sulphate	-	E5	5	% w/v	<i>myo</i> -Inositol	-
A6	2	M	Lithium sulphate	Partial	E6	10	% w/v	<i>myo</i> -Inositol	-
A7	20	% v/v	Ethylene glycol	-	E7	6	% w/v	L-Rhamnose monohydrate	-
A8	40	% v/v	Ethylene glycol	Yes	E8	12	% w/v	L-Rhamnose monohydrate	-
A9	80	% v/v	Ethylene glycol	Yes	E9	25	% w/v	L-Rhamnose monohydrate	-
A10	15	% v/v	1,2-Propandiol	-	E10	20	% v/v	Tetraethylene glycol (TTEG)	-
A11	30	% v/v	1,2-Propandiol	Partial	E11	40	% v/v	Tetraethylene glycol (TTEG)	Partial
A12	60	% v/v	1,2-Propandiol	Yes	E12	80	% v/v	Tetraethylene glycol (TTEG)	Yes
B1	15	% v/v	Glycerol	-	F1	20	% v/v	PEG 200	-
B2	30	% v/v	Glycerol	Partial	F2	40	% v/v	PEG 200	Yes
B3	60	% v/v	Glycerol	Yes	F3	80	% v/v	PEG 200	Yes
B4	15	% v/v	1,5-Pentanediol	-	F4	20	% v/v	Pentaethylene glycol	-
B5	30	% v/v	1,5-Pentanediol	Partial	F5	40	% v/v	Pentaethylene glycol	Yes
B6	60	% v/v	1,5-Pentanediol	Yes	F6	80	% v/v	Pentaethylene glycol	Yes
B7	20	% v/v	1,2,4-Butanetriol	-	F7	20	% v/v	D-Sorbitol	-
B8	40	% v/v	1,2,4-Butanetriol	Yes	F8	40	% v/v	D-Sorbitol	-
B9	80	% v/v	1,2,4-Butanetriol	Yes	F9	60	% v/v	D-Sorbitol	-
B10	20	% v/v	1,2-Hexanediol	-	F10	20	% v/v	Hexaethylene glycol	-
B11	40	% v/v	1,2-Hexanediol	-	F11	40	% v/v	Hexaethylene glycol	-
B12	80	% v/v	1,2-Hexanediol	-	F12	80	% v/v	Hexaethylene glycol	Yes
C1	20	% w/v	1,6-Hexanediol	-	G1	20	% v/v	PEG 300	-
C2	40	% w/v	1,6-Hexanediol	Partial	G2	40	% v/v	PEG 300	Partial
C3	80	% w/v	1,6-Hexanediol	Yes	G3	80	% v/v	PEG 300	Yes
C4	20	% v/v	MPD	-	G4	15	% w/v	Sucrose	-
C5	40	% v/v	MPD	Yes	G5	30	% w/v	Sucrose	-
C6	80	% v/v	MPD	Yes	G6	60	% w/v	Sucrose	Yes
C7	12	% w/v	<i>meso</i> -Erythritol	-	G7	20	% w/v	Maltitol	-
C8	25	% w/v	<i>meso</i> -Erythritol	-	G8	40	% w/v	Maltitol	-
C9	50	% w/v	<i>meso</i> -Erythritol	Yes	G9	80	% w/v	Maltitol	Yes
C10	20	% w/v	Trimethylolpropane	-	G10	12	% w/v	D-Trehalose	-
C11	40	% w/v	Trimethylolpropane	Yes	G11	25	% w/v	D-Trehalose	-
C12	80	% w/v	Trimethylolpropane	Yes	G12	50	% w/v	D-Trehalose	Yes
D1	20	% w/v	1,2,6-Hexanetriol	-	H1	20	% v/v	PEG 400	-
D2	40	% w/v	1,2,6-Hexanetriol	Yes	H2	40	% v/v	PEG 400	Partial
D3	80	% w/v	1,2,6-Hexanetriol	Yes	H3	80	% v/v	PEG 400	Yes
D4	20	% w/v	D-Lyxose	-	H4	20	% v/v	Pentaerythritol propoxylate (5/4 PO/OH)	-
D5	40	% w/v	D-Lyxose	-	H5	40	% v/v	Pentaerythritol propoxylate (5/4 PO/OH)	Partial
D6	80	% w/v	D-Lyxose	Yes	H6	80	% v/v	Pentaerythritol propoxylate (5/4 PO/OH)	Yes
D7	10	% w/v	Xylitol	-	H7	15	% v/v	PEG 500 MME	-
D8	20	% w/v	Xylitol	-	H8	30	% v/v	PEG 500 MME	-
D9	40	% w/v	Xylitol	Yes	H9	60	% v/v	PEG 500 MME	Yes
D10	20	% w/v	D-Fructose	-	H10	15	% v/v	PEG 600	-
D11	40	% w/v	D-Fructose	Yes	H11	30	% v/v	PEG 600	Partial
D12	80	% w/v	D-Fructose	Yes	H12	60	% v/v	PEG 600	Yes

Unit 5A, R-Evolution @
The Advanced Manufacturing Park
Selden Way, Rotherham, S60 5XA



**Molecular
Dimensions**

enquiries@moleculardimensions.com
+44 (0)114 224 2257

calibre
SCIENTIFIC