

Pi-PEG Screen

Cat.-No.: CS-128

Number	Precipitant 1	Buffer	pH	Additive
1	22.5 % v/v PEG 200	50 mM Acetate	4.8	None
2	20.0 % v/v PEG 300	50 mM Acetate	5.2	None
3	20.0 % v/v PEG 350 MME	50 mM MES	5.6	None
4	20.0 % v/v PEG 400	50 mM MES	6.0	None
5	20.0 % v/v PEG 550 MME	50 mM ADA	6.4	None
6	20.0 % v/v PEG 600	50 mM ADA	6.8	None
7	17.5 % w/v PEG 1000	50 mM HEPES	7.1	None
8	17.5 % w/v PEG 1500	50 mM HEPES	7.3	None
9	15.0 % w/v PEG 2000	50 mM Tris	7.6	None
10	15.0 % w/v PEG 2000 MME	50 mM Tris	8.0	None
11	15.0 % w/v PEG 3000	50 mM Bicine	8.4	None
12	15.0 % w/v PEG 4000	50 mM Bicine	8.8	None
13	15.0 % w/v PEG 1500	50 mM Acetate	4.8	3.6 % w/v PEG 4000
14	12.9 % w/v PEG 2000	50 mM Acetate	5.2	6.4 % v/v PEG 200
15	12.9 % w/v PEG 2000 MME	50 mM MES	5.6	5.7 % v/v PEG 300
16	12.9 % w/v PEG 3000	50 mM MES	6.0	5.7 % v/v PEG 350 MME
17	12.9 % w/v PEG 4000	50 mM ADA	6.4	5.7 % v/v PEG 400
18	19.3 % v/v PEG 200	50 mM ADA	6.8	5.7 % v/v PEG 550 MME
19	17.1 % v/v PEG 300	50 mM HEPES	7.1	5.7 % v/v PEG 600
20	17.1 % v/v PEG 350 MME	50 mM HEPES	7.3	5.0 % w/v PEG 1000
21	17.1 % v/v PEG 400	50 mM Tris	7.6	5.0 % w/v PEG 1500
22	17.1 % v/v PEG 550 MME	50 mM Tris	8.0	4.3 % w/v PEG 2000
23	17.1 % v/v PEG 600	50 mM Bicine	8.4	4.3 % w/v PEG 2000 MME
24	15.0 % w/v PEG 1000	50 mM Bicine	8.8	4.3 % w/v PEG 3000

Pi-PEG Screen

Cat.-No.: CS-128

Number	Precipitant 1	Buffer	pH	Additive
25	14.3 % v/v PEG 350 MME	50 mM Acetate	4.8	11.4 % v/v PEG 300
26	14.3 % v/v PEG 400	50 mM Acetate	5.2	11.4 % v/v PEG 350 MME
27	14.3 % v/v PEG 550 MME	50 mM MES	5.6	11.4 % v/v PEG 400
28	14.3 % v/v PEG 600	50 mM MES	6.0	11.4 % v/v PEG 550 MME
29	12.5 % w/v PEG 1000	50 mM ADA	6.4	11.4 % v/v PEG 600
30	12.5 % w/v PEG 1500	50 mM ADA	6.8	10.0 % w/v PEG 1000
31	10.7 % w/v PEG 2000	50 mM HEPES	7.1	10.0 % w/v PEG 1500
32	10.7 % w/v PEG 2000 MME	50 mM HEPES	7.3	8.6 % w/v PEG 2000
33	10.7 % w/v PEG 3000	50 mM Tris	7.6	8.6 % w/v PEG 2000 MME
34	10.7 % w/v PEG 4000	50 mM Tris	8.0	8.6 % w/v PEG 3000
35	16.1 % v/v PEG 200	50 mM Bicine	8.4	7.1 % w/v PEG 4000
36	14.3 % v/v PEG 300	50 mM Bicine	8.8	12.9 % v/v PEG 200
37	8.6 % w/v PEG 2000 MME	50 mM Acetate	4.8	17.1 % v/v PEG 400
38	8.6 % w/v PEG 3000	50 mM Acetate	5.2	17.1 % v/v PEG 550 MME
39	8.6 % w/v PEG 4000	50 mM MES	5.6	17.1 % v/v PEG 600
40	12.9 % v/v PEG 200	50 mM MES	6.0	15.0 % w/v PEG 1000
41	11.4 % v/v PEG 300	50 mM ADA	6.4	15.0 % w/v PEG 1500
42	11.4 % v/v PEG 350 MME	50 mM ADA	6.8	12.9 % w/v PEG 2000
43	11.4 % v/v PEG 400	50 mM HEPES	7.1	12.9 % w/v PEG 2000 MME
44	11.4 % v/v PEG 550 MME	50 mM HEPES	7.3	12.9 % w/v PEG 3000
45	11.4 % v/v PEG 600	50 mM Tris	7.6	10.7 % w/v PEG 4000
46	10.0 % w/v PEG 1000	50 mM Tris	8.0	19.3 % v/v PEG 200
47	10.0 % w/v PEG 1500	50 mM Bicine	8.4	17.1 % v/v PEG 300
48	8.6 % w/v PEG 2000	50 mM Bicine	8.8	17.1 % v/v PEG 350 MME

Pi-PEG Screen

Cat.-No.: CS-128

Number	Precipitant 1	Buffer	pH	Additive 1
49	8.6 % v/v PEG 550 MME	50 mM Acetate	4.8	22.9 % v/v PEG 600
50	8.6 % v/v PEG 600	50 mM Acetate	5.2	20.0 % w/v PEG 1000
51	7.5 % w/v PEG 1000	50 mM MES	5.6	20.0 % w/v PEG 1500
52	7.5 % w/v PEG 1500	50 mM MES	6.0	17.1 % w/v PEG 2000
53	6.4 % w/v PEG 2000	50 mM ADA	6.4	17.1 % w/v PEG 2000 MME
54	6.4 % w/v PEG 2000 MME	50 mM ADA	6.8	17.1 % w/v PEG 3000
55	6.4 % w/v PEG 3000	50 mM HEPES	7.1	14.3 % w/v PEG 4000
56	6.4 % w/v PEG 4000	50 mM HEPES	7.3	25.7 % v/v PEG 200
57	9.6 % v/v PEG 200	50 mM Tris	7.6	22.9 % v/v PEG 300
58	8.6 % v/v PEG 300	50 mM Tris	8.0	22.9 % v/v PEG 350 MME
59	8.6 % v/v PEG 350 MME	50 mM Bicine	8.4	22.9 % v/v PEG 400
60	8.6 % v/v PEG 400	50 mM Bicine	8.8	22.9 % v/v PEG 550 MME
61	4.3 % w/v PEG 4000	50 mM Acetate	4.8	25.0 % w/v PEG 1500
62	6.4 % v/v PEG 200	50 mM Acetate	5.2	21.4 % w/v PEG 2000
63	5.7 % v/v PEG 300	50 mM MES	5.6	21.4 % w/v PEG 2000 MME
64	5.7 % v/v PEG 350 MME	50 mM MES	6.0	21.4 % w/v PEG 3000
65	5.7 % v/v PEG 400	50 mM ADA	6.4	17.9 % w/v PEG 4000
66	5.7 % v/v PEG 550 MME	50 mM ADA	6.8	32.1 % v/v PEG 200
67	5.7 % v/v PEG 600	50 mM HEPES	7.1	28.6 % v/v PEG 300
68	5.0 % w/v PEG 1000	50 mM HEPES	7.3	28.6 % v/v PEG 350 MME
69	5.0 % w/v PEG 1500	50 mM Tris	7.6	28.6 % v/v PEG 400
70	4.3 % w/v PEG 2000	50 mM Tris	8.0	28.6 % v/v PEG 550 MME
71	4.3 % w/v PEG 2000 MME	50 mM Bicine	8.4	28.6 % v/v PEG 600
72	4.3 % w/v PEG 3000	50 mM Bicine	8.8	25.0 % w/v PEG 1000

Pi-PEG Screen

Cat.-No.: CS-128

Number	Precipitant 1	Buffer	pH	Additive 1
73	2.5 % w/v PEG 1000	50 mM Acetate	4.8	25.7 % w/v PEG 2000 MME
74	2.5 % w/v PEG 1500	50 mM Acetate	5.2	25.7 % w/v PEG 3000
75	2.1 % w/v PEG 2000	50 mM MES	5.6	21.4 % w/v PEG 4000
76	2.1 % w/v PEG 2000 MME	50 mM MES	6.0	38.6 % v/v PEG 200
77	2.1 % w/v PEG 3000	50 mM ADA	6.4	34.3 % v/v PEG 300
78	2.1 % w/v PEG 4000	50 mM ADA	6.8	34.3 % v/v PEG 350 MME
79	3.2 % v/v PEG 200	50 mM HEPES	7.1	34.3 % v/v PEG 400
80	2.9 % v/v PEG 300	50 mM HEPES	7.3	34.3 % v/v PEG 550 MME
81	2.9 % v/v PEG 350 MME	50 mM Tris	7.6	34.3 % v/v PEG 600
82	2.9 % v/v PEG 400	50 mM Tris	8.0	30.0 % w/v PEG 1000
83	2.9 % v/v PEG 550 MME	50 mM Bicine	8.4	30.0 % w/v PEG 1500
84	2.9 % v/v PEG 600	50 mM Bicine	8.8	25.7 % w/v PEG 2000
85	None	50 mM Acetate	4.8	25.0 % w/v PEG 4000
86	None	50 mM Acetate	5.2	45.0 % v/v PEG 200
87	None	50 mM MES	5.6	40.0 % v/v PEG 300
88	None	50 mM MES	6.0	40.0 % v/v PEG 350 MME
89	None	50 mM ADA	6.4	40.0 % v/v PEG 400
90	None	50 mM ADA	6.8	40.0 % v/v PEG 550 MME
91	None	50 mM HEPES	7.1	40.0 % v/v PEG 600
92	None	50 mM HEPES	7.3	35.0 % w/v PEG 1000
93	None	50 mM Tris	7.6	35.0 % w/v PEG 1500
94	None	50 mM Tris	8.0	30.0 % w/v PEG 2000
95	None	50 mM Bicine	8.4	30.0 % w/v PEG 2000 MME
96	None	50 mM Bicine	8.8	30.0 % w/v PEG 3000