


# Ligand fitting puzzles 2

Judit Debreczeni




Ligands meeting

07 12 2016



# Problem 1.




- mixed ligand population issue workflow
- smi → pdb, cif with grade, acedrg and pyrogen
  - acedrg fell over on the smi file if it contained description (does the newer version fix this?)
- basic ligand fitting in coot
  - coot misplaced the ligand – not enough conformers tested?
  - ligand expert is now default – good!
  - manual intervention: fix chirality – not sure where it went wrong!
- extra density which can't be water due to close proximity – probably missing atom..



# Problem 1.





• extra density which can't be water due to close proximity – probably missing atom...  
PDB 1LEE and 1LF2: CRYSTAL STRUCTURE OF PLASMEPSIN FROM *P. FALCIPARUM* IN COMPLEX WITH INHIBITOR RS367

- useful to know what's in the tube
- 2 similar ligands R36 and R37 – probably have a mixture of both in 1LEE.
- model this as two different ligands at occ=0.5
- generate cif dictionary for second form (did this by modifying the smi file and re-running pyrogen/acedrg), different (non-default) three letter code
- need to merge two cif files for remlac to run (buster accepts two cif files separately which is easier for the user)
- Lidia falls over!



## Problem 2.

- high resolution, straight-forward ligand fitting
- smi → pdb, cif with grade, acedrg and pyrogen
- ligand fitting in coot – worked well except Me-piperazine – wrong conformation here, needed manual intervention.
- piperazine is also probably protonated – close proximity to acid side chain, but this is missed in Lidia even though we do have a 2.9Å contact. It is also quite easy to fit this group the wrong way round.



## Confidentiality Notice

This file is private and may contain confidential and proprietary information. If you have received this file in error, please notify us and remove it from your system and note that you must not copy, distribute or take any action in reliance on it. Any unauthorized use or disclosure of the contents of this file is not permitted and may be unlawful. AstraZeneca PLC, 1 Francis Crick Avenue, Cambridge Biomedical Campus, Cambridge, CB2 0AA, UK, T: +44(0)203 749 5000, [www.astrazeneca.com](http://www.astrazeneca.com)

